



User's Quick Install Guide

LifeSentry

Instructions for Installing Your New Personal Emergency Response System

Your Kit Includes:



Rechargeable AA Batteries



Lithium Ion Pendant Batteries



Telephone Cord



AC Adapter



Lanyard



Belt Clip



Wrist Band

Equipment Reference:

Base Unit Bottom View



Master On/Off Switch

Battery Backup Door

Base Unit Top View



Power Light indicates current power condition:
OFF: No power to unit
ON SOLID: Unit has Power
ON FLASHING: Running on battery power

Pendant Battery Charger Door

Info Light alerts you to fault conditions:
OFF: Normal Operation
FLASHING: Problem with one or more pendants

Base Unit Back View



Line In Connects to telephone jack

To House Phone

Programming Switch This should be set to Normal Operation "Emergency Call Mode"

Power Cord Insert plug and twist right to lock

Pendant

Front View



Emergency Button

Back View



Battery Test, System Test and Call Disconnect

Contact Information:

To reach an Advanced Product Support Representative who can help with step-by-step installation please call 843-369-7803.

User's Quick Install Guide

Step 1: Place Batteries in Unit



(A) Place the spare pendant Lithium Ion battery in the charger on the front of the Base Unit. (Make sure charging light turns RED.)



(B) Place the Green Rechargeable batteries in the compartment on the bottom of the Base Unit. Make sure to turn ON/OFF switch to ON. You will hear the announcement "Running on Battery Power." Unit will do so until the AC Adapter power supply is plugged in.



(C) Place one Lithium Ion battery into the pendant by unscrewing the cap and inserting FLAT SIDE UP. As soon as cap is tightened, the red light flashes as the pendant does a self-check and announces the battery status.

Step 2: Plug in the Telephone and Power Cords



(A) Plug the AC Adaptor power cord into the power jack and twist right to lock.



(B) Plug one end of the supplied phone cord into the "Line-in" jack. Connect the other end into the wall (or modem) phone jack.



(C) Plug the telephone into the "To House Phone" jack.

Step 3 (REQUIRED): Contact HTC to confirm set up is complete

Contact HTC at 843-369-2222 to confirm that your unit has been set up and to review your contact list. Please have your unit number ready when you call. It can be found on the bottom of the Base Unit.

Step 4: Test the unit



Press the blue emergency button for four seconds to send a test alarm. The green light will flash. You will hear dialing and other sounds as the alarm is sent. When the monitoring station answers, just tell them you are testing. Once this is done set-up will be complete.

Attention: We recommend that you test the battery regularly by pressing the gray battery test button on the back of the Pendant. Also, be sure to keep your spare Pendant battery charged. Simply place the spare Pendant battery in the charger and keep it there for when you need to change batteries. Don't worry, it is ok to leave the battery in the charger until you need it. The charger is designed to keep the battery from overcharging. The battery in the Pendant will stand by for many months. However, we recommend that you exchange batteries monthly.